

BIO

Insight into health

5000L

Color Doppler Ultrasound System

BIO 5000L

Color Doppler Ultrasound System



Wuxi Biomedical Technology Co., Ltd.
Address: Room 40B, No.801 Hongqiao road,
Wuxi, Jiangsu, China 214072
Tel/Fax: 86-510-85572789
Email: sales@cnbiomedical.com





BIO 5000L

Color Doppler Ultrasound System



The portable Color Doppler system A5 with superior image, high performance and low cost is your right choice. Ricso makes efforts to let each person can get high-quality medical services. A5 is the one that can meet your various requirement with an affordable price. User-friendly interface can make doctors work more efficiently.


Lightweight Design



12 inch LED high-resolution display



Main machine weight 5.5kg



Thickness of mainframe only 13.8cm

Full range of configuration

- | A variety of transducers
- | Built-in lithium battery
- | Professional multi-purpose cart, visiting package
- | HDMI and multi-USB interface, data transmission is more convenient



Micro-convex probe



Convex probe





Linear probe



Transvaginal probe



AQ-Sonic, E-Crystal

① Select different options depending on the different characteristics of different tissues. Speckle and noise reduction can improve the imaging quality.



② Three simultaneous displays: 2D, Color doppler and PW is simultaneously displayed in real time, in favor of real-time observation and comparison.

③ Auto trace: calculating the real-time data automatically, making diagnosis faster and more convenient.



④ Ultra-wideband beamformer.

⑤ Intelligent imaging optimization technology.

