

BI+ BM6+

Portable Digital Microscope Dermatoscope

Features:


- n 3 inch TFT color display, rotated 360°
- n 8 LEDs with brightness adjustable
- n 5M sensor & professional microscopic lens
- n A metal stand fixes the dermatoscope
- n 200x magnification(up to 500x by digital)
- n Video, photo and sound capture
- n Works with Li-ion rechargeable battery or power supply.
- n Micro SD card optional, download via USB.
- n Measurement by software (1/1000mm)





Product specifications:

Image sensor	5M pixels
Lens	High quality microscopic lens
Display	3 inches TFT panel, 4:3 ratio
Magnification	20x-200x, up to 500x by digital
Still Image resolution	12M (4030x3024), 9M, 5M, 3M, 1.3M, VGA
Video capture resolution	VGA, QVGA
Frame rate	30fps@VGA
Light source	8 LEDs with brightness adjustable
Power source	(1) Li-ion 3.7V/1050mAh rechargeable battery (2) Adaptor: Input: 100-240V; Output: 5V 1A (3) USB
TV output	To any monitor with TV-in
Storage	Support Micro-SD card up to 32GB SDHC (not included)
OSD Language	English/Germany/French/Spanish/Italian/ Portuguese/Japanese/Chinese
Measurement	By Micro Capture Pro Software (Window XP/Vista/7/8, Mac 10.7 and above)
Software Language	English/Germany/French/Spanish
Dimension	130x103x29mm
Weight	195g(With battery)

Pictures



Oberse the mole



Oberse hair

Wuxi Biomedical Technology Co., Ltd.

www.portable-ultrasound-scanner.com

Add:Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel:86-510-85572789

Fax:86-510-85572789

Email:sales@cnbiomedical.com

Mobilephone:86-13382885910

BI@ is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.